

## EMPTY GLASS VIAL INSPECTION



Free Glass Fragment Inside Detected
Resolution Until 150 micron
Inspection of Inside and FULL Surface
Dimensional Analisys
Class A Designed
Full Automatic Motorized Changeover
• Easy To Use

**EGVI-** Performs the complete control of the Full Surface of Vials.

It assure the quality of the empty vials before the filling operation.

Bottom, Mouth, And Sidewall Lateral Surface are inspected

It detect free glass fragments inside the vials, detect foreign particles inside or outside vials, detect critical glass defect, perform precise dimensional analysis.


EGVI can be installed directly in Sterile Environments, before the filling machine.

Easy to install thanks to linear designed.


Full Automatic and motorized changeover.

Modular design, small footprint


Specifically designed for sterile room installation.



Layout Bottom & Mouth Control



Layout Bottom , Mouth and Lateral SideWall Control



PRODUCT	Glass Vials
CHANGEOVER	Full Automatic
GLASS TYPE	Tube or molded
SPEED	Up 36.000 Vials /min
RESOLUTION	70 micron

